Langen, den 10.06.2022

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Bibliographie

Quellenangaben

1. Deerain et al Dec 21 2022
Assessment of the Analytical Sensitivity of 10 Lateral Flow Devices against the SARS-CoV-2 Omicron Variant. J Clin Microbiol. 2022 Feb 16;60(2)
<https://journals.asm.org/doi/10.1128/jcm.02479-21>

⁹ https://ec.europa.eu/health/system/files/2022-05/covid-19_rat_common-list_en.pdf

¹⁰ <https://eur-lex.europa.eu/eli/reg/2017/746/oj>

¹¹ https://lexparency.de/eu/32017R0746/ART_100/

2. Apurv Soni et al Febr 2022
Comparison of Rapid Antigen Tests' Performance between Delta (B.1.61.7; AY.X) and Omicron (B.1.1.529; BA1) Variants of SARS-CoV-2: Secondary Analysis from a Serial Home Self-Testing Study
www.medrxiv.org/content/10.1101/2022.02.27.22271090v1
3. James Regan et al Dec 2021
Detection of the omicron variant virus with the Abbott BinaxNow SARS-CoV-2 Rapid Antigen Assay
www.medrxiv.org/content/10.1101/2021.12.22.21268219v1
4. Statens Serum Institut, Copenhagen, DK. 08 Mar 2022
Testing of SARS-CoV-2 rapid antigen tests detection of variants (Delta and Omicron BA.1 and BA.2)
<https://en.ssi.dk/-/media/arkiv/subsites/covid19/diagnostik/afprvning-af-sars-cov-2-antigentests-for-pvisning-af-varianter.pdf?la=en>
5. Goderski G et al (RIVM, NL) 16 Dec 2021
Technical evaluation of SARS-CoV-2 antigen self-tests with Omicron variant
www.rivm.nl/sites/default/files/2021-12/Technical-evaluation-of-SARS-CoV-2-Self-test-with-omicron-variant_Final.pdf
6. UK Health Security Agency. Technical briefing 32. 17 Dec 2021
SARS-CoV-2 variants of concern and variants under investigation in England – Technical briefing 32
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/1042046/Technical_Briefing_32.pdf
7. Zachary A. Weishampel et al 9 Feb 2022
OraSure IntelliSwab™ Rapid Antigen Test Performance with the SARS-CoV-2 Variants of Concern—Alpha, Beta, Gamma, Delta, and Omicron
www.mdpi.com/1999-4915/14/3/543
8. John Schrom et al 10 Jan 2022
Direct Comparison of SARS Co-V-2 Nasal RT- PCR and Rapid Antigen Test (BinaxNOW™) at a Community Testing Site During an Omicron Surge
www.medrxiv.org/content/10.1101/2022.01.08.22268954v4
9. Molenkamp, R & Igloi, Z. 2022
Evaluation of Antigen rapid test and PCR test to Omicron variant
www.erasmusmc.nl/-/media/erasmusmc/pdf/1-themaspecifiek/viroscience/2021-evaluation-omicron-in-pcr-and-ag-assays.pdf

10. Barbara L Goodall et al. 21 Jan 2022
Investigating sensitivity of nasal or throat (ISNOT): A combination of both swabs increases sensitivity of SARS-CoV-2 rapid antigen tests
www.medrxiv.org/content/10.1101/2022.01.18.22269426v1
11. Kanjilal et al. 19 Feb 2022
Performance of three rapid antigen tests against the SARS-CoV-2 Omicron variant
www.medrxiv.org/content/10.1101/2022.02.17.22271142v1
12. de Michelena P et al 22 Feb 2022
Real-life performance of a COVID-19 rapid antigen detection test targeting the SARS-CoV-2 nucleoprotein for diagnosis of COVID-19 due to the Omicron variant.
[www.journalofinfection.com/article/S0163-4453\(22\)00110-4/fulltext](https://www.journalofinfection.com/article/S0163-4453(22)00110-4/fulltext)
13. Bayart, J.-L. et al 01 Mar 2022
Analytical Sensitivity of Six SARS-CoV-2 Rapid Antigen Tests for Omicron Versus Delta Variant
www.preprints.org/manuscript/202203.0010/v1
14. Eckerle et al 17 Jan 2022
Sensitivity of SARS-CoV-2 antigen-detecting rapid tests for Omicron variant
www.medrxiv.org/content/10.1101/2021.12.18.21268018v2
15. Cochrane COVID-19 Diagnostic Test Accuracy Group, 24 Mar 2021
Rapid, point-of-care antigen and molecular-based tests for diagnosis of SARS-CoV-2 infection
www.cochranelibrary.com/cdsr/doi/10.1002/14651858.CD013705.pub2/full
16. Denkinger et al Aug 2021 eCollection
Accuracy of novel antigen rapid diagnostics for SARS-CoV-2: A living systematic review and meta-analysis
<https://pubmed.ncbi.nlm.nih.gov/34383750/>
17. Navarro et al; 22 Feb 2022
Real-life performance of a COVID-19 rapid antigen detection test targeting the SARS-CoV-2 nucleoprotein for diagnosis of COVID-19 due to the Omicron variant; The Journal of Infection 84(5): e64-e66
<https://europepmc.org/article/pmc/pmc8863328>
18. Related Articles zu Nr. 18
[www.journalofinfection.com/article/S0163-4453\(22\)00110-4/fulltext#relatedArticles](https://www.journalofinfection.com/article/S0163-4453(22)00110-4/fulltext#relatedArticles)
19. Premsrirut et al; 05 Jan 2022
Discordant SARS-CoV-2 PCR and Rapid Antigen Test Results When

Infections: A December 2021 Occupational Case Series

www.medrxiv.org/content/10.1101/2022.01.04.22268770v1

20. Klimek-Tulwin Rapid Reviews 03 Mar 2022

Review2: Sensitivity of SARS-CoV-2 antigen-detecting rapid tests for Omicron variant

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[1] Deerain et al; Feb 16 2022

Assessment of the Analytical Sensitivity of 10 Lateral Flow Devices against the SARS-CoV-2 Omicron Variant. J Clin Microbiol; 60(2)

<https://journals.asm.org/doi/10.1128/jcm.02479-21>

The analytical sensitivities of the 10 antigen kits were similar for both the Delta and Omicron variants.

[2] Apurv Soni et al Feb 2022

Comparison of Rapid Antigen Tests' Performance between Delta (B.1.61.7; AY.X) and Omicron (B.1.1.529; BA1) Variants of SARS-CoV-2: Secondary Analysis from a Serial Home Self-Testing Study

www.medrxiv.org/content/10.1101/2022.02.27.22271090v1

The performance of Ag-RDT is not inferior among individuals infected with the SARS-CoV-2 Omicron variant as compared to the Delta variant. The improvement in sensitivity of Ag-RDT noted with serial testing is consistent between Delta and Omicron variant.

[3] James Regan et al Dec 2021

Detection of the omicron variant virus with the Abbott BinaxNow SARS-CoV-2 Rapid Antigen Assay

www.medrxiv.org/content/10.1101/2021.12.22.21268219v1

All Omicron and delta specimens with concentrations of 100,000 copies/swab or greater were positive by the BinaxNow Assay, a concentration similar to previously reported limits of detection for this assay.

[4] Statens Serum Institut, Copenhagen, DK. 08.03.2022

Testing of SARS-CoV-2 rapid antigen tests detection of variants (Delta and Omicron BA.1 and BA.2)

<https://en.ssi.dk/-/media/arkiv/subsites/covid19/diagnostik/afprvning-af-sars-cov-2-antigentests-for-pvisning-af-varianter.pdf?la=en>

The testing of the 23 tested antigen tests shows that all tests can detect Delta and Omicron BA.1 and BA.2 on a comparable level with the wild-type (Wuhan).

[5] Goderski G et al (RIVM, NL) 16 Dec 2021

Technical evaluation of SARS-CoV-2 antigen self-tests with Omicron variant
www.rivm.nl/sites/default/files/2021-12/Technical-evaluation-of-SARS-CoV-2-Self-test-with-omicron-variant_Final.pdf

In conclusion, the performance of all [7] self-tests is similar and not affected by the Omicron variant, and therefore meeting the pre-set criteria.

[6] UK Health Security Agency. Technical briefing 32. 17 Dec 2021

SARS-CoV-2 variants of concern and variants under investigation in England – Technical briefing 32

https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/1042046/Technical_Briefing_32.pdf

Five different Ag RDT Lateral Flow Tests detect Omicron as accurately as Delta.

[7] Zachary A. Weishampel et al 9 Febr 2022

OraSure IntelliSwab™ Rapid Antigen Test Performance with the SARS-CoV-2 Variants of Concern—Alpha, Beta, Gamma, Delta, and Omicron

www.mdpi.com/1999-4915/14/3/543

Ultimately, the OraSure IntelliSwab™ COVID-19 Rapid Test showed no decrease in sensitivity between the ancestral SARS-CoV-2 strain and any VOCs including Omicron.

[8] John Schrom et al 10 Jan 2022

Direct Comparison of SARS Co-V-2 Nasal RT- PCR and Rapid Antigen Test (BinaxNOW™) at a Community Testing Site During an Omicron Surge

<https://www.medrxiv.org/content/10.1101/2022.01.08.22268954v4>

731 people found that the Abbott BinaxNOW rapid tests performed about as well with Omicron as they did with other variants when people were symptomatic and had high viral load

[9] Molenkamp, R & Igloi, Z. 2022

Evaluation of Antigen rapid test and PCR test to Omicron variant

www.erasmusmc.nl/-/media/erasmusmc/pdf/1-themaspecifiek/viroscience/2021-evaluation-Omicron-in-pcr-and-ag-assays.pdf

Three rapid Ag-tests showed comparable sensitivity Omicron to Delta

[10] Barbara L Goodall et al. 21 Jan 2022

Investigating sensitivity of nasal or throat (ISNOT): A combination of both swabs increases sensitivity of SARS-CoV-2 rapid antigen tests

www.medrxiv.org/content/10.1101/2022.01.18.22269426v1

Omicron detection, comparable to previous evaluations using other SARS-CoV-2 lineages

[11] Kanjilal et al 19 Febr 2022

Performance of three rapid antigen tests against the SARS-CoV-2 Omicron variant

<https://www.medrxiv.org/content/10.1101/2022.02.17.22271142v1>

We found no change in the analytic sensitivity of all three RADTs for detection of Omicron versus Delta, but noted differences in performance between assays

[12] de Michelena P et al 22 Febr 2022

Real-life performance of a COVID-19 rapid antigen detection test targeting the SARS-CoV-2 nucleoprotein for diagnosis of COVID-19 due to the Omicron variant.

[www.journalofinfection.com/article/S0163-4453\(22\)00110-4/fulltext](https://www.journalofinfection.com/article/S0163-4453(22)00110-4/fulltext)

In our experience, the clinical performance of the Panbio™ COVID-19 assay for Omicron variant was comparable to previous reports for the Wuhan-Hu-1 G614 variant (100% specificity and sensitivity of 81.4% in non-vaccinated adult patients with a clinical course of <5 days)

[13] Bayart, J.-L. et al 01 Mar 2022

Analytical Sensitivity of Six SARS-CoV-2 Rapid Antigen Tests for Omicron Versus Delta Variant

www.preprints.org/manuscript/202203.0010/v1

Sensitivity across a wide range of Ct values (13.5 to 35.7; median = 21.3) were comparable and ranged from 70.0% to 77.1% for Delta strains and from 69.6% to 78.3% for Omicron strains