



**FOR NCA REVIEW ONLY –
DO NOT PUBLISH**

BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins – Switzerland
Tél: +41 21 556 30 Fax:
+41 21 556 30 99
www.BD.com

XX September 2022

URGENT: FIELD SAFETY NOTICE – IDS-22-4562

BD Sabouraud Agar with Penicillin and Streptomycin

REF: 254451 Lot Numbers: See Table 1

Type of Action: Product Removal

**Attention: Clinical Personnel, Risk Managers, Laboratory Manager,
Laboratory Technician**

This letter contains important information which requires your **immediate** attention.

Dear Customer,

BD is conducting a Field Safety Corrective Action to remove specific lots of BD Sabouraud Agar with Penicillin and Streptomycin as listed in Table 1 below. According to our distribution records, your organisation may have received the impacted product.

Product Code (REF)	Product Description	Lot Number	Expiry Date
254451	Sabouraud agar with Penicillin and Streptomycin (20 plates)	2152635	04/10/2022
		2153805	29/09/2022
		2193936	10/11/2022

Table 1: Impacted product

Description of the problem

As a result of internal testing, BD has identified a potential supplier issue with one of the plate's impregnated antibiotics, streptomycin, which may impact the effectivity of the agent to inhibit bacterial growth over the shelf life of the plate.

This product removal is limited to the product code / lot numbers listed in Table 1 above. No other product codes or lot numbers are affected.



**FOR NCA REVIEW ONLY –
DO NOT PUBLISH**

BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins – Switzerland
Tél: +41 21 556 30 Fax:
+41 21 556 30 99
www.BD.com

Clinical risk

This defect could result in the erroneous identification of pathogens in clinical specimens. These potential outcomes are due to the possibility of inappropriate bacterial growth within the plate, which can be incorrectly identified as a fungal organism or can inhibit the growth and detection of fungal organisms. In these instances, erroneous results may lead to a delay in diagnosis, treatment of a patient with an antimicrobial that is inappropriate or ineffective, or clinical worsening of disease state.

To date there has been no adverse events worldwide related to this issue.

Actions taken by BD

BD is working with the supplier to identify the root cause and implement corrective actions to prevent recurrence of this issue.

Actions for Customers:

1. Inspect your inventory, locate and quarantine any plates of the impacted lot numbers. Make a note of these lot numbers, then destroy all affected product. For plates destroyed, credit will be provided.
2. Circulate this Field Safety Notice to all those within your organisation that may use the BD Sabouraud Agar with Penicillin and Streptomycin as listed in the Table 1 above.
3. If you have further distributed the product, please identify those facilities, notify them at once and ensure they follow the instructions listed in this notification.
4. There are no additional follow-up activities recommended for patient samples processed with these plates.
5. If you experience any issues with the BD Sabouraud Agar with Penicillin and Streptomycin, please report as a complaint as per your normal process.
6. Complete and sign the Customer Response Form on page 4 and return it to **<<insert contact details here>>** as soon as possible or no later than **21st October 2022**, indicating the following:
 - the number of plates destroyed, **OR**
 - that your organisation does not have any impacted plates left in inventory.

Note: If you no longer have or use the product, it is still important that you return the Customer Response Form for our reconciliation purposes.



**FOR NCA REVIEW ONLY –
DO NOT PUBLISH**

BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins – Switzerland
Tél: +41 21 556 30 Fax:
+41 21 556 30 99
www.BD.com

Contact reference person

If you have any questions about this, please contact your local BD representative or the local BD office on <<insert telephone details here>> or e-mail <<insert contact email address here>>.

We confirm that the appropriate regulatory agencies have been informed of these actions.

BD is committed to advancing the world of health. Our primary objectives are patient safety and user safety and providing you with quality products. We apologise for the inconvenience this situation may cause you and thank you in advance for helping BD to resolve this matter as quickly and effectively as possible.

Sincerely,

A handwritten signature in black ink, appearing to read 'L. Darrock'.

Lorna Darrock
Sr. Manager, Post Market Quality
EMEA Quality



**FOR NCA REVIEW ONLY –
DO NOT PUBLISH**

BD Switzerland Sàrl
Terre Bonne Park – A4
Route de Crassier 17
1262 Eysins – Switzerland
Tél: +41 21 556 30 Fax:
+41 21 556 30 99
www.BD.com

Customer Response Form - IDS-22-4562
BD Sabouraud Agar with Penicillin and Streptomycin

REF: 254451 Lot Numbers: See Table 1

Return to << insert email address >> as soon as possible or **no later than the 21st October 2022.**

I confirm this notice has been read, understood and that all recommended actions have been implemented as required.

Tick the appropriate box below

We do not have any of the affected product as listed in Table 1 in our possession.

OR

We have the following plates of the affected product as listed in Table 1 in our possession and I confirm that the plates have been destroyed (*Please complete the table below with the number of plates destroyed for each lot. Credit will only be provided on completion and return of this form.*)

REF:	Lot Number/s:	Plates destroyed <i>(insert unit quantity below)</i>
254451	2152635	
	2153805	
	2193936	

Account/Organisation Name:	
Department <i>(if applicable):</i>	
Address:	
Postcode:	City:
Contact Name:	Job Title:
Contact Telephone Number:	Contact E-mail Address:
Name of your supplier for this product <i>(if not direct from BD*):</i>	
Signature:	Date:

This form must be returned to BD before this action can be considered closed for your account.

**If you were forwarded this Field Safety Notice via a distributor/3rd party, please return your completed form to that organisation for reconciliation purposes.*